

Goat Anti-Rabbit IgG Heavy and Light Chain Antibody, Fluorescein Conjugated

ORDERING INFORMATION

Catalog No: 81301

Size: 1ml at 1mg/ml

SPECIFICATION SUMMARY

Antigen: Rabbit IgG

Host Species: Goat

Stabilizers: 0.2% Bovine Serum Albumin (BSA)

Preservatives: 0.09% Sodium Azide

Buffer: Phosphate-buffered saline (PBS)

Purification: Immunoaffinity chromatography

Conjugation: F/P ratio is 4.85

SPECIFICITY

By immunoelectrophoresis and ELISA this antibody reacts specifically with rabbit IgG and with light chains common to other rabbit immunoglobulins. Reactivity with non-immunoglobulin serum proteins is not detected. This antibody may cross-react with IgG from other species.

APPLICATIONS

Most Applications: 1:200 – 1:2,000

STORAGE AND STABILITY

Store at 2-8°C. Shelf life is one (1) year from date of receipt.

For in vitro research use only. Not for use in diagnostic or therapeutic procedures.