

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Goat Anti-Mouse IgG Fc Fragment Specific Antibody, FITC Conjugated

ORDERING INFORMATION

Catalog No: 83301

Size: 1ml at 1mg/ml

SPECIFICATION SUMMARY

Antigen: Mouse IgG Fc fragment

Host Species: Goat

Stabilizers: 0.2% Bovine Serum Albumin (BSA)

Preservatives: 0.09% sodium azide

Buffer: Phosphate-buffered saline (PBS)

Purification: Immunoaffinity chromatography

Conjugation: F/P ratio = 4.4

SPECIFICITY

By immunoelectrophoresis and ELISA this antibody reacts specifically with mouse IgG-Fc. Cross-reactivity with IgA and IgM is negligible. This antibody may cross-react with IgG from other species.

APPLICATIONS

Immunohistochemistry: 1:200-1:2,000

Immunocytochemistry: 1:200-1:2,000

Prepare working dilutions immediately before use. These are recommended dilutions.

User should determine optimal dilutions for their applications.

STORAGE AND STABILITY

Store at 2-8°C. Shelf life is one (1) year from date of receipt.

For in vitro research use only. Not for use in diagnostic or therapeutic procedures.