

Goat anti Rabbit IgG, Fc Fragment Specific Antibody, Peroxidase Conjugated

ORDERING INFORMATION

Catalog No: 84101

Size: 1ml at 1mg/ml

SPECIFICATION SUMMARY

Antigen: Rabbit IgG Fc fragment

Host Species: Goat

Stabilizers: 0.2% Bovine Serum Albumin

Preservatives: 0.1% Pro-Clean 400

Buffer: Phosphate-buffered saline

Purification: Immunoaffinity chromatography

Molar enzyme:antibody ratio: 4:1

SPECIFICITY

By immunoelectrophoresis and ELISA this antibody reacts specifically with rabbit IgG Fc fragment. This antibody may cross-react with IgG of other species.

APPLICATIONS

ELISA: 1:10,000 – 1:100,000

Western blot: 1:1,000 – 1:30,000

Immunohistochemistry: 1:200 - 1:500

Immunocytochemistry: 1:200 - 1:500

These are recommended concentrations. Enduser should determine optimal concentrations for their applications.

DILUTION INSTRUCTIONS

Dilute in PBS or in medium that is identical to that used in the assay system.

STORAGE AND STABILITY

Store at 2-8°C. Expires 1 year from date of receipt.

For in vitro investigational use only. Not for use in diagnostic or therapeutic applications.