

## Goat Anti-Human IgE Antibody

### **ORDERING INFORMATION**

**Catalog No:** 87091

**Size:** 1ml at 1mg/ml

### **SPECIFICATION SUMMARY**

**Antigen:** Purified human IgE

**Host Species:** Goat

**Stabilizers:** None

**Preservatives:** 0.1% sodium azide

**Purification:** Immunoaffinity chromatography

### **SPECIFICITY**

By immunoelectrophoresis and ELISA this antibody reacts specifically with human IgE.

### **APPLICATIONS**

ELISA Coating:	1:100-1:500
Primary:	1:1,000-1:30,000
Western Blot Colorimetric	1:1,000-1:10,000
Chemiluminescent:	1:1,000-1:30,000
Immunohistochemistry:	1:200-1:2,000
Immuocytochemistry:	1:200-1:2,000

These are recommended dilutions only. Enduser should determine optimal dilutions for their applications.

### **STORAGE AND STABILITY**

Store at 2-8°C. Expires 2 years from date of receipt.

**For *in vitro* investigational use only. Not for therapeutic or diagnostic procedures.**