

## Goat anti Human IgG, Heavy and Light Chain Specific Antibody, FITC Conjugated

### **ORDERING INFORMATION**

**Catalog No:** 87305

**Size:** 2.0 mg, lyophilized powder.

### **SPECIFICATION SUMMARY**

**Antigen:** Human IgG

**Host Species:** Goat

**Buffer:** 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6

**Stabilizers:** 15 mg/ml Bovine Serum Albumin

**Preservatives:** 0.05% Sodium Azide

**Purification:** Immunoaffinity chromatography

### **SPECIFICITY**

The antibody reacts with the heavy chains and light chains associated with human immunoglobulin G.

### **MOST APPLICATIONS**

Begin diluting at: 1:50 – 1:200

### **DILUTION INSTRUCTIONS**

Reconstitute with 1.5 mL distilled water, or a combination of 0.75 ml of distilled water and 0.75 ml of glycerol if stored frozen.

### **STORAGE AND STABILITY**

Store at 2-8°C until reconstituted. Upon reconstitution, this product should be stored frozen in 50% glycerol. Expires 1 year from the date of reconstitution.

**For *in vitro* research use only.**