

QED Bioscience Inc.

ADVANCED RESEARCH TECHNOLOGIES

Goat Anti-Human IgE Antibody, Biotin Conjugated

ORDERING INFORMATION

Catalog No: 87491

Size: 1.0ml, 1mg/ml

SPECIFICATION SUMMARY

Antigen: Purified human IgE

Host Species: Goat

Stabilizers: 0.2% bovine serum albumin

Preservatives: 0.09% sodium azide

Purification: Immunoaffinity chromatography

SPECIFICITY

By immunoelectrophoresis and ELISA, this antibody reacts specifically with human IgE.

APPLICATIONS

ELISA 1: 10,000-1:100,000

Immunoblotting 1: 10,000-1:100,000

These are recommended dilutions. Enduser should determine optimal dilution for their application.

DILUTION INSTRUCTIONS

Dilute in PBS or medium identical to that used in assay system. Prepare working dilutions immediately before use.

STORAGE AND STABILITY

Store at 2-8°C. stable for 1 year from date of receipt.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.