

ChromaLink™ One-Shot Antibody Biotinylation Kit

Catalog Number: **QB-9007-009K** Size: **Kit** Price: **\$150.00**

- **Complete Kit:** All necessary reagents included to biotinylate 100ug of antibody
- **Built-in Quantitative Signal:** Incorporated biotin quantified by simple absorbance – no HABA or fluoro-reporter assay needed – only a simple UV scan needed to measure
- **Safe and Economical:** Water soluble, mild on antibodies – no DMF needed



No HABA or
Fluoro-Reporter Assay




UV-traceable antibody
biotinylation



Foolproof protocol
on a flash drive

ChromaLink Biotin One-Shot Antibody Labeling Kit offers the complete solution to the quantitative biotinylation of 100ug of antibody . The kit incorporates traceable biotin and allows for direct spectrophotometric verification of biotin incorporation. A simple Å280 and Å354 measurement determines the number of biotins per antibody molecule.



Product	Catalog#	Description	Size/Quantity
ChromaLink™ Biotin One-Shot Antibody Labeling Kit	QB-9007-009	Traceable biotinylation	For 100ug of antibody
		Sulfo-ChromaLink™ Biotin	8.0ug/1
		1X Modification Buffer pH 7.4	1.5ml/1
		1X PBS pH 7.2	1.5ml/1
		Collection tubes	1.5ml/4
		Zeba™ Spin Columns	0.5ml/2
		Biotinylated IgG Control	100ug/1
		1M Tris HCl pH 8.9	100ul/1
		0.5GB Flash drive	1