

Anti-Monkey IgG Monoclonal Antibodies

ORDERING INFORMATION

Catalog No.: 14801 (clone 5C12.D4)
and 14804 (clone 4D8.B10)

Size: 1 ml ascites or 1 mg Protein G-
purified antibody in PBS, pH 7.4

SPECIFICATION SUMMARY

Antigen: Chromatographically purified
Rhesus monkey IgG.

Host Species: Mouse

Antibody Class: IgG1

Stabilizers: None

Preservatives: None. Available on
request.

ACTIVITY

	Ig mg/ml*	ELISA Titer**
14801	25.0	>1 : 102,400
14804	4.0	1: 200

* Determined by radial immunodiffusion
for several lots of ascites.

**ELISA is an indirect assay against
monkey IgG on the solid phase. Titers
represent endpoint dilutions that achieve
absorbance values ≥ 0.1 .

SPECIFICITY

These antibodies recognize monkey IgG
and do not cross-react with human,
bovine, equine or ovine IgG.

APPLICATIONS

ELISA

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is
identical to that used in the assay
system.

STORAGE AND STABILITY

These antibodies are stable for at least
one (1) year at -20°C to -70°C . Store
product in appropriate aliquots to avoid
multiple freeze-thaw cycles.

*For in vitro investigational use only. Not
for use in therapeutic or diagnostic
procedures.*