

Anti-TRICK2 (DR5, TRAIL-R2) Antibody

ORDERING INFORMATION

Catalog No.: 11091

Size: 100ug IgG in PBS, pH 7.4, purified by immunoaffinity chromatography.

BACKGROUND

Apoptosis is induced by cytokines in the TNF family through their death domain containing receptors. TRICK2 (TRAIL receptor inducer of cell killing 2) is a member of the TNF receptor family and contains the conserved death domain. TRICK2 is expressed at high levels in lymphocytes and spleen cells and in many types of tumor cells. Overexpression of TRICK2 rapidly induces cell death by apoptosis. TRICK2 has been shown to function as a receptor for TRAIL and can inhibit TRAIL-induced killing of lymphocytic cells. Other names for TRICK2 include DR5, Apo2, TRAIL-R2, and KILLER.

SPECIFICATION SUMMARY

Immunogen: Synthetic peptide corresponding to amino acids 409-427 (LGERLAKQKIEDHLLSSGK, isoform 1) of human TRICK2 precursor conjugated to KLH at the N-terminus. The same sequence is found at amino acids 380-398 in isoform 2.

Host Species: Rabbit

Stabilizers: None

Preservatives: None. Available on request.

SPECIFICITY

This antibody recognizes human TRICK2.

APPLICATIONS

Immunoblotting: use a 1:500 dilution. A band of approximately 57kDa is detected.

Positive control: whole cell lysate from HL60, HeLa, or Jurkat cells.

These are recommended concentrations. Enduser should determine optimal concentrations for their applications.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C. Avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.