

ATM Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 13012

Size: 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

SPECIFICATION SUMMARY

Antigen: Synthetic peptide representing a portion of the protein encoded in part by exons 37 and 38 (LocusLink ID 472).

Host Species: Goat

Stabilizers: None

Preservatives: 0.1% Sodium Azide

SPECIFICITY

This antibody reacts with human and mouse. Antibody is affinity purified.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.

APPLICATIONS

Western Blot: 0.2 - 1 ug/ml.

Immunoprecipitation: 10 – 20 ug/mg lysate

STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.