

ATM Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 13014

Size: 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

SPECIFICATION SUMMARY

Antigen: Synthetic peptide representing a portion of the protein encoded within exon 53 (LocusLink ID 472).

Host Species: Goat

Stabilizers: None

Preservatives: 0.1% Sodium Azide

SPECIFICITY

This antibody reacts with human and mouse. Antibody is affinity purified.

APPLICATIONS

Western Blot: 0.1 – 0.4 ug/ml.

Immunoprecipitation: 10 – 20 ug/mg lysate

STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.