

ATR Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 13016

Size: 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

SPECIFICATION SUMMARY

Antigen: Synthetic peptide representing a portion of the protein encoded within exon 5 (LocusLink ID 545).

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.1% Sodium Azide

SPECIFICITY

This antibody reacts with human and mouse. Antibody is affinity purified.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.

APPLICATIONS

Western Blot: 0.033 – 0.067 ug/ml.

STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.