

## CKII alpha Polyclonal Antibody

### **ORDERING INFORMATION**

**Catalog No.:** 13028

**Size:** 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

### **SPECIFICATION SUMMARY**

**Antigen:** Synthetic peptide representing a portion of the protein located within exons 12 and 13 (LocusLink ID 1457).

**Host Species:** Rabbit

**Stabilizers:** None

**Preservatives:** 0.1% Sodium Azide

### **SPECIFICITY**

This antibody reacts with human. Other species have not been tested. Antibody is affinity purified.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*

### **APPLICATIONS**

*Western Blot:* 0.3 – 1 ug/ml.

*Immunoprecipitation:* Use Catalog Number 13030.

### **STORAGE AND STABILITY**

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.