

Anti-Rabbit Immunoglobulin Monoclonal Antibodies

ORDERING INFORMATION

Catalog No.	Clone No.	MAb Subtype	Size	Library Pack No.	100ug/clone
14601	RIG10-110.1	IgG1	100ug, 500ug	146101	2 clones
14602	RIG10-120.22	IgG1	100ug, 500ug		

Format: Protein G-purified antibody in PBS, pH 7.4.

SPECIFICATION SUMMARY

Antigen: Pool of purified rabbit immunoglobulin from normal sera.

Host Species: Mouse

Specificity: #14601 recognizes rabbit IgG and IgM; #14602 recognizes rabbit IgG. These antibodies do not cross-react with immunoglobulin of the following species: bovine, canine, chicken, donkey, equine, feline, rat, sheep, porcine, guinea pig, goat, or human.

APPLICATIONS

These antibodies have been qualified for use in ELISA to detect rabbit immunoglobulin.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

These antibodies are stable for at least one (1) year at -20°C to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.