

Ketanserin Monoclonal Antibody

ORDERING INFORMATION

Catalog No.: 17101

Size: 500 ug Protein G-purified antibody in PBS, pH 7.4.

SPECIFICATION SUMMARY

Antigen: Ketanserin conjugated via carboxymethylxime to keyhole limpet hemocyanin (KLH).

Host Species: Mouse

Antibody Class: IgG1

Stabilizers: None

Preservatives: None. Available on request.

ACTIVITY AND SPECIFICITY

This antibody has a Kd of 2 nM for ketanserin.

APPLICATIONS

This antibody may be used in immunoassays for detection of ketanserin. Other applications are under investigation.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

These antibodies are stable for at least one (1) year at -20°C to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.

For in vitro investigational use only.

Not for use in therapeutic or diagnostic procedures.