

Hepatitis C Virus Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 18604

Size: 100 ug affinity-purified antibody
in PBS, pH 7.4

SPECIFICATION SUMMARY

Antigen: Peptide corresponding to
amino acids 9-21 of Hepatitis C Virus
(HCV) core protein C1.

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.02% sodium azide.

SPECIFICITY

This antibody recognizes the core
protein C1 of Hepatitis C Virus.

APPLICATIONS

Immunoblotting: use at 1:500-1:1,000.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is
identical to that used in the assay
system.

STORAGE AND STABILITY

This antibody is stable for at least one
(1) year at -70°C. Avoid multiple freeze-
thaw cycles.

*For in vitro investigational use only. Not
for use in therapeutic or diagnostic
procedures.*