

Hepatitis C Virus Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 18605

Size: 100 ug affinity-purified antibody
in PBS, pH 7.4

SPECIFICATION SUMMARY

Antigen: Peptide corresponding to
amino acids 33-43 of Hepatitis C Virus
(HCV) core protein C3.

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.02% sodium azide.

SPECIFICITY

This antibody reacts with the 44 kD core
protein C3, the envelope protein, and
E2/NS1 of Hepatitis C Virus.

APPLICATIONS

Immunoblotting: use at 1:500-1:1,000.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is
identical to that used in the assay
system.

STORAGE AND STABILITY

This antibody is stable for at least one
(1) year at -70°C. Avoid multiple freeze-
thaw cycles.

*For in vitro investigational use only. Not
for use in therapeutic or diagnostic
procedures.*