

## Human Plasminogen Monoclonal Antibody

### **ORDERING INFORMATION**

**Catalog No.:** 21401-21403

**Size:** 500 ug Protein A-purified antibody in PBS, pH 7.4.

### **SPECIFICATION SUMMARY**

**Antigen:** Purified human plasminogen.

**Host Species:** Mouse

**Antibody Class:** 21401 = IgG2b;  
21402 = IgG1; 21403 = IgG2a.

**Stabilizers:** None

**Preservatives:** 0.1% sodium azide.

### **SPECIFICITY**

These antibodies are specific for human plasminogen.

### **APPLICATIONS**

These antibodies may be used in immunoassays to detect/quantitate human plasminogen. They may also be used to affinity purify plasminogen and in immunohistochemistry on formalin-fixed/paraffin embedded tissues.

### **DILUTION INSTRUCTIONS**

Dilute in PBS or medium which is identical to that used in the assay system.

### **STORAGE AND STABILITY**

These antibodies are stable for at least one (1) year at -20°C to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*