


**recomLine EBV IgG [Avidity] [IgA]**  
**recomLine CMV IgM**

**New Antigens**  
**ZEBRA, BZLF1**

Strip immunoassay with recombinant antigens for detection of IgG, IgA, and IgM antibodies to Epstein-Barr virus (EBV).

Epstein-Barr virus (EBV), an ubiquitously occurring Herpes virus, can cause symptoms of infectious mononucleosis on primary infection. As a result of lifelong persistence of this virus, reactivations can occur, especially in immunocompromised individuals.

*recomLine* EBV utilizes recombinant EBV antigens on nitrocellulose strips for detection of IgG, IgA, and IgM antibodies to multiple EBV antigens on one strip.



**Product #4572** *recomLine* EBV IgG [Avidity] [IgA] (20 tests per kit)

**Product #4573** *recomLine* EBV IgM (20 tests per kit)

For more information visit <http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.  
 Not intended for diagnostic procedures.*