

**recomWell Yersinia IgG**  
**recomWell Yersinia IgM**


ELISA format with recombinant antigens for the detection of IgG and IgM antibodies to *Yersinia enterocolitica* and *Yersinia pseudotuberculosis*.

*Yersinia enterocolitica* has become increasingly important as an intestinal bacterium that is pathogenic for humans. This organism is transmitted orally through contaminated water or food. Typical symptoms of an acute *Yersinia* infection are diarrhea, stomach ache, and fever.

*recomWell Yersinia* assays use recombinant *Yersinia* outer membrane proteins (YOPs) of only *Yersinia* strains that are pathogenic for humans.

**Product #4604** *recomWell Yersinia IgG*  
(96 tests per kit)

**Product #4605** *recomWell Yersinia IgM*  
(96 tests per kit)



For more information visit  
<http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.  
Not intended for diagnostic procedures.*