

## MIKROGEN DIAGNOSTIK





ELISA format with recombinant antigens for detection of IgG and IgM antibodies to Hepatitis E virus (HEV).

Hepatitis E virus (HEV) is one of the leading causes of fecal-oral acquired hepatitis worldwide, usually resulting from contaminated drinking water. In developed countries, HEV is usually transmitted through infected pork that was not sufficiently cooked. HEV can also be transmitted through blood products, blood transfusions, and organ donation. HEV infections often remain asymptomatic, but acute HEV infections present with symptoms similar to those of Hepatitis A.

*recom*Well HEV IgG and *recom*Well IgM utilize recombinant antigens of the HEV genotypes 1 and 3.

Product #5004 recomWell HEV IgG (96 tests per kit)

Product #5005 recomWell HEV IgM (96 tests per kit)



For more information visit <u>http://www.mikrogen.de/english/deutschland</u>/home.html

Available in the U.S. For Research Use Only. Not intended for diagnostic procedures.