

Goat Anti Chicken IgY, Heavy and Light Chain Specific Antibody, FITC Conjugated

ORDERING INFORMATION

Catalog No: 86301

Size: 1ml at 1.0mg/ml

SPECIFICATION SUMMARY

Antigen: Chicken IgY

Host Species: Goat

Stabilizers: 0.2% BSA

Preservatives: 0.09% sodium azide

Buffer: Phosphate buffered saline

Purification: Immunoaffinity chromatography

Conjugation: F/P ratio is 3 to 7.

SPECIFICITY

By immunoelectrophoresis and ELISA this antibody reacts specifically with chicken IgY and with light chains common to other chicken immunoglobulins. Reactivity with other non-immunoglobulin serum proteins is not detected. This antibody may cross-react with IgG of other species.

APPLICATIONS

Immunohistochemistry: 1:200-1:2,000

Immunocytochemistry: 1:200-1:2,000

STORAGE AND STABILITY

Store at 2-8°C. Expires 1 year from date of receipt.

For *in vitro* investigational use only. Not for diagnostic or therapeutic procedures.