

## CKII alpha Polyclonal Antibody

### ORDERING INFORMATION

**Catalog No.:** 13030

**Size:** 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

### SPECIFICATION SUMMARY

**Antigen:** Synthetic peptide representing a portion of the protein located within exon 14 (LocusLink ID 1457).

**Host Species:** Rabbit

**Stabilizers:** None

**Preservatives:** 0.1% Sodium Azide

### SPECIFICITY

This antibody reacts with human. Other species have not been tested. Antibody is affinity purified.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*

### APPLICATIONS

*Western Blot:* Use Catalog Number 13026 or 13028.

*Immunoprecipitation:* 4 –10 ug/mg lysate.

### STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.