

CKII alpha' Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 13032

Size: 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

SPECIFICATION SUMMARY

Antigen: Synthetic peptide representing a portion of the protein located within exon 11 (LocusLink ID 1459).

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.1% Sodium Azide

SPECIFICITY

This antibody reacts with human. Other species have not been tested. Antibody is affinity purified.

APPLICATIONS

Western Blot: 0.2 – 1 ug/ml.

Immunoprecipitation: 4–10 ug/mg lysate.

STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.