

p70S6 Kinase Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 13038

Size: 100 ug IgG in 0.1 ml Tris-citrate/phosphate, pH 7-8

SPECIFICATION SUMMARY

Antigen: Synthetic peptide representing a portion of the C-terminus of the protein (LocusLink ID 6197).

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.1% Sodium Azide

SPECIFICITY

This antibody reacts with rat. Other species have not been tested; mouse and human p70S6 are identical to the site of the immunizing sequence. Antibody is affinity purified.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.

APPLICATIONS

Western Blot: 0.05-0.4 ug/ml.

STORAGE AND STABILITY

Store at 2 - 8°C. Antibody is stable at 2 - 8°C for 1 year.