

## Anti-Canine Immunoglobulin Monoclonal Antibodies

### ORDERING INFORMATION

Catalog No.	Clone No.	MAB Subtype	Size	Library Pack No.	100ug/clone
14001	DIG11-124.1	IgG1	100ug, 500ug	140101	All 4 clones
14002	DIG11-127.2	IgG1	100ug, 500ug		
14003	DIG12-223.3	IgG1	100ug, 500ug		
14004	DIG12-255.5	IgG1	100ug, 500ug		

**Format:** Protein G-purified antibody in PBS, pH 7.4.

### SPECIFICATION SUMMARY

**Antigen:** Pool of purified canine immunoglobulin from normal sera.

**Host Species:** Mouse

**Specificity:** These antibodies recognize canine immunoglobulin and do not cross-react with immunoglobulin of the following species: bovine, equine, feline, rabbit, sheep, porcine, donkey, guinea pig, goat, or human.

### APPLICATIONS

These antibodies have been qualified for use in ELISA to detect canine immunoglobulin.

### DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

### STORAGE AND STABILITY

These antibodies are stable for at least one (1) year at -20°C to -70°C.

Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*