

Anti-Bovine Immunoglobulin Monoclonal Antibodies

ORDERING INFORMATION

Catalog No.	Clone No.	MAB Subtype	Size	Library Pack No.	100ug/clone
14101	BIG10-101.6	IgG1	100ug, 500ug	141101	All 4 clones
14102	BIG10-108.5	IgG1	100ug, 500ug		
14103	BIG10-123.1	IgG1	100ug, 500ug		
14104	BIG10-137.4	IgG1	100ug, 500ug		

Format: Protein G-purified antibody in PBS, pH 7.4.

SPECIFICATION SUMMARY

Antigen: Pool of purified bovine immunoglobulin from normal sera.

Host Species: Mouse

Specificity: These antibodies recognize bovine immunoglobulin and do not cross-react with immunoglobulin of the following species: canine, equine, feline, rabbit, sheep, porcine, donkey, guinea pig, goat, or human.

APPLICATIONS

These antibodies have been qualified for use in ELISA to detect bovine immunoglobulin.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

These antibodies are stable for at least one (1) year at -20°C to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.