

***Citrobacter koseri* Polyclonal Antibody**

ORDERING INFORMATION

Catalog No.: 18730

Size: 100 ul

Format: Antiserum

SPECIFICATION SUMMARY

Antigen: *Citrobacter koseri*

Host Species: Rabbit

Preservatives: None

Affinity: 0.2×10^{-3} g/l

SPECIFICITY

Antibody specificity determined by a competitive ELISA.

APPLICATIONS

EIA, ImmunoHisto/Cytochemistry:

1:15,000 – 1:30,000

STORAGE AND STABILITY

Stable for 2 years when stored at -20°C.

Stable for 2 months at 4°C when stored with 0.1% merthiolate or a 50/50 mixture of PBS and glycerol.

Compound	Cross-Reactivity Ratios
<i>Hafnia alvei</i>	$>10^{-2}$ g/l
<i>Morganella morganii</i>	$>10^{-2}$ g/l
<i>Pseudomonas aeruginosa</i>	$>10^{-2}$ g/l
<i>Pseudomonas putida</i>	$>10^{-2}$ g/l
<i>Proteus mirabilis</i>	$>10^{-2}$ g/l
<i>Klebsiella pneumoniae</i>	$>10^{-2}$ g/l

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.