

## Full-Length Retinoblastoma Protein (p110<sup>RB</sup>)

### **ORDERING INFORMATION**

**Catalog No.:** 3108

**Size:** 35ug

**Concentration:** 1.3 mg/ml.

### **STORAGE AND STABILITY**

This product is stable for at least six (6) months at -70°C or below.

*For in vitro investigational use only.*

*Not for use in diagnostic or therapeutic procedures.*

### **SPECIFICATION SUMMARY**

**Source:** Genetically engineered *E. coli*.

**Molecular Weight:** 110,000 daltons.

**Purity:** >90%.

**Diluent:** 20 mM sodium phosphate

200 mM sodium chloride

1 mM EDTA

10% glycerol

pH 7.5

Sterile-filtered

**Preservatives:** None.

**Biologically Active:** Yes.

**Activity Assay:** [<sup>3</sup>H]-Thymidine incorporation anti-proliferation bioassay.