

## Retinoblastoma (Rb) phosphoserine 612 Polyclonal Antibody

### ORDERING INFORMATION

**Catalog No.:** 3113

**Size:** 25ug

**Format:** Sequential immunoaffinity-purified antibody in phosphate buffer, pH 7.4, 0.05% sodium azide, 1mg/ml BSA.

### SPECIFICATION SUMMARY

**Antigen:** Phosphopeptide region of human Rb protein containing serine 612.

**Host Species:** Rabbit

**Specificity:** This antibody recognizes human Rb protein phosphorylated at serine 612. Reactivity with other species has not been investigated.

### APPLICATIONS

**Immunoblotting:** use at 0.15-0.5ug/ml. Endusers should determine optimal concentration for their application.

### STORAGE AND STABILITY

Store at -80°C. Prepare appropriate aliquots to avoid multiple freeze-thaw cycles.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*