

Retinoblastoma (Rb) phosphoserine 811 Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 3116

Size: 25ug

Format: Sequential immunoaffinity-purified antibody in phosphate buffer, pH 7.4, 0.05% sodium azide, 1mg/ml BSA.

SPECIFICATION SUMMARY

Antigen: Phosphopeptide region of human Rb protein containing serine 811.

Host Species: Rabbit

Specificity: This antibody recognizes human Rb protein phosphorylated at serine 811. Reactivity with other species has not been investigated.

APPLICATIONS

Immunoblotting: use at 0.50-1.5ug/ml. Endusers should determine optimal concentration for their application.

STORAGE AND STABILITY

Store at -80°C. Prepare appropriate aliquots to avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.