

BAD (Phospho-Ser136) Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 43068

Format: 100ul at 1.0mg/ml in PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Affinity-purified on phosphopeptide; non-phosphopeptide-reactive antibodies were removed by chromatography on non-phosphorylated peptide.

BACKGROUND

The protein encoded by the *BAD* gene is a member of the Bcl-2 family whose members are regulators of programmed cell death. BAD protein positively regulates cell apoptosis by forming heterodimers with Bcl-xL and Bcl-2 and reversing their death repressor activity. Proapoptotic activity of this protein is regulated through its phosphorylation. Protein kinases AKT and MAP kinase, as well as protein phosphatase calcineurin, were found to be involved in the regulation of this protein. Alternative splicing of this gene results in two transcript variants which encode the same isoform.

SPECIFICATION SUMMARY

Antigen: Peptide sequence that includes phosphorylation site of serine 136 (S-R-S(p)-A-P) derived from mouse BAD and conjugated to KLH.


Host Species: Rabbit

Specificity: This antibody detects endogenous human, mouse, and rat BAD only when phosphorylated at serine 136.

Accession no.: Q61337, NP_031548.1


APPLICATIONS

Immunoblotting: use at dilution of 1:500-1:1,000. A band of ~23kDa is detected.



Detection of BAD (phospho-Ser136) in extracts of HT29 cells untreated or treated with UV.

Immunohistochemistry: use at dilution of 1:50-1:100.



Detection of BAD (phospho-Ser136) in paraffin-embedded human breast carcinoma tissue.

These are recommended working dilutions. Enduser should determine optimal dilutions for their application.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C. Can be stored at 4°C for short-term.

For in vitro investigational use only. Not intended for therapeutic or diagnostic applications.