

recomWell Yersinia IgG
recomWell Yersinia IgA
recomWell Yersinia IgM

ELISA format with recombinant antigens for the detection of IgG, IgA, and IgM antibodies to *Yersinia enterocolitica* and *Yersinia pseudotuberculosis*.


Yersinia enterocolitica has become increasingly important as an intestinal bacterium that is pathogenic for humans. This organism is transmitted orally through contaminated water or food. Typical symptoms of an acute *Yersinia* infection are diarrhea, stomach ache, and fever.

recomWell Yersinia assays use recombinant *Yersinia* outer membrane proteins (YOPs) of only *Yersinia* strains that are pathogenic for humans.

Product #4604 *recomWell Yersinia IgG*
(96 tests per kit)

Product #4605 *recomWell Yersinia IgM*
(96 tests per kit)

Product #4606 *recomWell Yersinia IgA*
(96 tests per kit)



For more information visit
<http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.
Not intended for diagnostic procedures.*