

recomLine Yersinia IgG 2.0
recomLine Yersinia IgA [IgM] 2.0

NEW


Strip-Immunoassay with recombinant antigens for the detection of IgG, IgM, and IgA antibodies to *Yersinia enterocolitica* and *Yersinia pseudotuberculosis*.

Yersinia enterocolitica has become increasingly important as an intestinal bacterium that is pathogenic for humans. This organism is transmitted orally through contaminated water or food. Typical symptoms of an acute *Yersinia* infection are diarrhea, stomach ache, and fever.

recomLine Yersinia assays use recombinant *Yersinia* outer membrane proteins (YOPs) and adhesins of only *Yersinia* strains that are pathogenic for humans.

Product #4672 recomLine Yersinia IgG 2.0
(20 tests per kit)

Product #4673 recomLine Yersinia IgA [IgM] 2.0
(20 tests per kit)



For more information visit

<http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.
 Not intended for diagnostic procedures.*