


**recomWell Helicobacter IgG**  
**recomWell Helicobacter IgA**

ELISA format with recombinant antigens for detection of IgG and IgA antibodies to *Helicobacter pylori*.

The gram-negative, spiral-shaped bacterium *Helicobacter pylori* was described by Warren and Marshall in 1983. Since then, the connection between chronic Type B gastritis and *H. pylori* has been well-documented. In fact, there is a causal relationship between *H. pylori* and almost all duodenal and stomach ulcers. Moreover, *H. pylori* is associated with an increased risk of gastric adenocarcinomas. Apart from strains of different virulence, host factors may also be responsible for this phenomenon.

*H. pylori* infections induce IgG and IgA antibody responses; usually only low levels of IgM antibodies are found.



**Product #4704 recomWell Helicobacter IgG**  
(96 tests per kit)

**Product #4705 recomWell Helicobacter IgA**  
(96 tests per kit)

For more information visit

<http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.  
Not intended for diagnostic procedures.*