


recomWell Helicobacter IgG
recomWell Helicobacter IgA

ELISA format with recombinant antigens for detection of IgG and IgA antibodies to *Helicobacter pylori*.

The gram-negative, spiral-shaped bacterium *Helicobacter pylori* was described by Warren and Marshall in 1983. Since then, the connection between chronic Type B gastritis and *H. pylori* has been well-documented. In fact, there is a causal relationship between *H. pylori* and almost all duodenal and stomach ulcers. Moreover, *H. pylori* is associated with an increased risk of gastric adenocarcinomas. Apart from strains of different virulence, host factors may also be responsible for this phenomenon.

H. pylori infections induce IgG and IgA antibody responses; usually only low levels of IgM antibodies are found.



Product #4704 recomWell Helicobacter IgG
(96 tests per kit)

Product #4705 recomWell Helicobacter IgA
(96 tests per kit)

For more information visit

<http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.
Not intended for diagnostic procedures.*