



**NEW** *recomWell* HEV IgG  
*recomWell* HEV IgM

ELISA format with recombinant antigens for detection of IgG and IgM antibodies to Hepatitis E virus (HEV).

Hepatitis E virus (HEV) is one of the leading causes of fecal-oral acquired hepatitis worldwide, usually resulting from contaminated drinking water. In developed countries, HEV is usually transmitted through infected pork that was not sufficiently cooked. HEV can also be transmitted through blood products, blood transfusions, and organ donation. HEV infections often remain asymptomatic, but acute HEV infections present with symptoms similar to those of Hepatitis A.

*recomWell* HEV IgG and *recomWell* IgM utilize recombinant antigens of the HEV genotypes 1 and 3.



**Product #5004** *recomWell* HEV IgG  
(96 tests per kit)

**Product #5005** *recomWell* HEV IgM  
(96 tests per kit)

For more information visit

<http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.  
Not intended for diagnostic procedures.*