


recomWell Treponema IgG
recomWell Treponema IgM

ELISA format with recombinant antigens for the detection of IgG and IgM antibodies to *Treponema pallidum*.

Treponema pallidum subsp. pallidum is the causative agent of syphilis, a chronic disease occurring only in humans. The disease progresses through stages. After a short propagation period at the site of infection (primary stage), spreading of the organism follows (secondary stage). After the decline of acute symptoms, *Treponema pallidum* can persist for years (latency stage) and then cause later complications (tertiary syphilis).

Recombinant *T. pallidum* antigens are used in the *recomWell* Treponema assays.



Product #5104 recomWell Treponema IgG
(96 tests per kit)

Product #5105 recomLineTreponema IgM
(96 tests per kit)

For more information visit
<http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.
Not intended for diagnostic procedures.*