


NEW *recomLine* CMV IgG [Avidity]
recomLine CMV IgM

Line immunoassay with recombinant antigens for detection of IgG and IgM antibodies to Cytomegalovirus (CMV) and to determine CMV-IgG avidity.

Infections with Cytomegalovirus (CMV) normally cause mild symptoms and can even remain asymptomatic. However, immunocompromised individuals and pregnant women are at high risk for suffering damage from a CMV infection; approximately 40% transmit the virus to the fetus.

recomLine CMV is a line immunoassay that utilizes recombinant phase-specific and avidity CMV antigens (MIKROGEN patented method) to detect IgG and IgM antibodies and the avidity of IgG antibodies.



Product #5572 *recomLine* CMV IgG [Avidity]
(20 tests per kit)

Product #5573 *recomLine* CMV IgM
(20 tests per kit)

For more information visit
<http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.
Not intended for diagnostic procedures.*