

DENTT Polyclonal Antibody

ORDERING INFORMATION

Catalog No.: 56316

Size: 100 ug affinity-purified antibody
in phosphate-buffered saline, pH 7.4.

SPECIFICATION SUMMARY

Antigen: Human DENTT protein.

Host Species: Rabbit

Stabilizers: None

Preservatives: 0.02% sodium azide

SPECIFICITY

This antibody recognizes human DENTT (Differentially Expressed Nucleolar TGF- β 1 Target) protein, a gene product expressed in several normal tissues with the highest level of expression in brain. DENTT is involved in responses mediated by TGF- β 1 and has distinct domains for direction to different subnuclear locations.

APPLICATIONS

Immunoblotting: use at 1:500-1:1,000 dilution. In immunoblots, a band of approximately 77 kD is detected.

DILUTION INSTRUCTIONS

Dilute in PBS or medium which is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least one (1) year at -20°C. Avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.