


recomWell Toxoplasma IgG **recomWell Toxoplasma IgM**

ELISA format with recombinant antigens for the detection of IgG and IgM antibodies to *Toxoplasma gondii*.

Toxoplasma gondii is an obligate intracellular, parasitic protozoan that causes toxoplasmosis. Found worldwide, *T. gondii* is capable of infecting virtually all warm-blooded animals, but felines such as domestic cats are the only known hosts in which the parasite can undergo sexual reproduction. Although mild, flu-like symptoms occasionally occur during the first few weeks following exposure, infection with *T. gondii* usually produces no readily observable symptoms in healthy human adults. However, in infants, HIV/AIDS patients, and others with weakened immunity, infection can cause serious and occasionally fatal toxoplasmosis.

Recombinant antigens of *T. gondii* tachyzoites and bradyzoites are used in the *recomWell* Toxoplasma assays.



Product #5904 *recomWell* Toxoplasma IgG
(96 tests per kit)

Product #5905 *recomWell* Toxoplasma IgM
(96 tests per kit)

For more information visit
<http://www.mikrogen.de/english/deutschland/home.html>

*Available in the U.S. For Research Use Only.
Not intended for diagnostic procedures.*