

recomLine Campylobacter IgG recomLine Campylobacter IgA

Increased
sensitivity
by use of p39


Strip immunoassay with recombinant antigens for detection of IgG or IgA antibodies to *Campylobacter jejuni* and *Campylobacter coli*

The genus *Campylobacter* comprises gram-negative, spiral-shaped, microaerophilic, mesophilic to thermophilic bacteria with bipolar flagella. Human *Campylobacter* infections are mainly food-associated intestinal infections with world-wide incidence. Infection by *Campylobacter jejuni* is much more frequent than *Campylobacter coli*.

recomLine Campylobacter kits detect IgG or IgA antibodies to specific recombinant antigens of *Campylobacter jejuni* and *Campylobacter coli*.

Product #6272 recomLine Campylobacter IgG
(20 tests per kit)

Product #6273 recomLine Campylobacter IgA
(20 tests per kit)



For more information visit

<http://www.mikrogen.de/english/deutschland/home.html>

Available in the U.S. For Research Use Only.
Not intended for diagnostic procedures.