

Goat Anti-Mouse IgA Antibody, HRP Conjugated

ORDERING INFORMATION

Catalog No.	Size
80405	1mg

Format: 1mg/ml in phosphate-buffered saline (PBS), 0.2% BSA, 0.1% Pro-Clean 400, conjugated to horseradish peroxidase. Molar enzyme:antibody ratio = 4:1.

SPECIFICATION SUMMARY

Antigen: Mouse IgA

Host Species: Goat

Specificity: Antibody is affinity-purified on mouse IgA. Cross-reactivity with mouse IgM, IgG1, IgG2a, IgG2b, IgG2c, IgG3, and IgE is not detectable.

APPLICATIONS

Immunoblotting: use at a dilution of 1:2,000-1:20,000.

Immunohistochemistry: use at a dilution of 1:200-1:2,000.

Immunocytochemistry: use at a dilution of 1:200-1:2,000.

ELISA: use at a dilution of 1:10,000-1:100,000.

These are recommended dilutions. Endusers should determine optimal dilutions for their applications. Prepare working dilutions immediately before use.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least 1 year at 2-8°C.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.