

Goat Anti-Mouse IgG1 Antibody, FITC Conjugated

ORDERING INFORMATION

Catalog No.	Size
83404	1.0mg

Format: 1mg/ml in phosphate-buffered saline (PBS), 0.2% BSA, 0.09% sodium azide, conjugated to fluorescein isothiocyanate. Fluorophore:antibody ratio = 4.9.

SPECIFICATION SUMMARY

Antigen: Mouse IgG1

Host Species: Goat

Specificity: Antibody is affinity-purified on mouse IgG1. Cross-reactivity with mouse IgA, IgM, IgG2a, IgG2b, IgG2c, IgG3, and IgE is not detectable.

APPLICATIONS

Immunohistochemistry: use at a dilution of 1:50-1:500.

Immunocytochemistry: use at a dilution of 1:50-1:500.

Immunofluorescence: use at a dilution of 1:50-1:500.

These are recommended dilutions. Endusers should determine optimal dilutions for their applications. Prepare working dilutions immediately before use.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least 1 year at 2-8°C.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.