

Goat Anti-Rabbit IgG-Fc Antibody, Biotin Conjugated

ORDERING INFORMATION

Catalog No: 84401

Size: 1ml at 1mg/ml

SPECIFICATION SUMMARY

Antigen: Rabbit IgG Fc fragment

Host Species: Goat

Stabilizers: 0.2% Bovine Serum Albumin (BSA)

Preservatives: 0.09% Sodium Azide

Buffer: Phosphate-buffered saline (PBS)

Purification: Immunoaffinity chromatography

SPECIFICITY

By immunoelectrophoresis and ELISA this antibody reacts specifically with rabbit IgG-Fc. Reactivity with rabbit Fab fragments is not detected. This antibody may cross-react with IgG from other species.

APPLICATIONS

ELISA: 1:10,000-1:100,000

Western Blot: 1:10,000-1:100,000

Prepare working dilutions immediately before use. These are recommended dilutions. User should determine optimal dilutions for their applications.

STORAGE AND STABILITY

Store at 2-8°C. Shelf life is one (1) year from date of receipt.

For in vitro research use only. Not for use in diagnostic or therapeutic procedures.