

Rabbit Anti-Chicken IgY Heavy and Light Chain Antibody, Peroxidase Conjugated

ORDERING INFORMATION

Catalog No.: 86403

Size: 1ml at 1mg/ml

SPECIFICATION SUMMARY

Antigen: Chicken IgY

Host Species: Rabbit

Format: Affinity-purified rabbit IgG in PBS, 0.2% BSA, 0.1% Pro-Clean 400

Molar enzyme/antibody ratio: 4:1

SPECIFICITY

By immunoblotting and ELISA this antibody reacts specifically with chicken IgY and with light chains common to all chicken immunoglobulins. Cross-reactivity with immunoglobulin of other species has not been evaluated.

APPLICATIONS

Immunoblotting: use at 1:1,000-1:30,000

Immunohistochemistry: use at 1:200-1:5,000

ELISA: use at 1:10,000-1:100,000

Prepare working dilutions immediately before use. These are recommended dilutions. Enduser should determine optimal dilutions for their applications.

STORAGE AND STABILITY

Store at 2-8°C. This antibody is stable for at least one (1) year from date of receipt.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.