

Goat anti Human IgM, Mu Chain Specific Antibody, Peroxidase Conjugated

ORDERING INFORMATION

Catalog No: 87120

Size: 2.0 ml, lyophilized powder.

SPECIFICATION SUMMARY

Antigen: Human IgM

Host Species: Goat

Buffer: 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6

Stabilizers: 15mg/ml Bovine Serum Albumin

Preservatives: None

Purification: Immunoaffinity chromatography

SPECIFICITY

The antibody reacts with the heavy chain associated with human immunoglobulin M.

APPLICATIONS

ELISA and Western Blot with Chromogenic Substrates: 1:5,000 – 1:100,000

Western Blot with ECL Substrates: 1:10,000 – 1:200,000

Immunohistochemistry: 1:500 – 1:5,000

DILUTION INSTRUCTIONS

Reconstitute with 2.0 mL distilled water, or a combination of 1.0 ml of distilled water and 1.0 ml of glycerol if stored frozen.

STORAGE AND STABILITY

Store at 2-8°C until reconstituted. Upon reconstitution, this product should be stored frozen in 50% glycerol. Expires 1 year from the date of reconstitution.

For *in vitro* research use only.