

Goat Anti-Rat IgG heavy and light chain antibody

ORDERING INFORMATION

Catalog No.	Size
88201	1.0mg

Format: 1mg/ml in phosphate-buffered saline (PBS), 0.09% sodium azide.

SPECIFICATION SUMMARY

Antigen: Rat IgG

Host Species: Goat

Specificity: Antibody is affinity-purified on rat IgG. This antibody may cross-react with IgG of other species.

APPLICATIONS

Immunoblotting: use at a dilution of 1:1,000-1:30,000

Immunohistochemistry: use at a dilution of 1:200-1:2,000.

Immunocytochemistry: use at a dilution of 1:200-1:2,000.

ELISA: for detection use at a dilution of 1:1,000-1:30,000, for coating use at a dilution of 1:100-1:500.

These are recommended dilutions. Endusers should determine optimal dilutions for their applications. Prepare working dilutions immediately before use.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least 1 year at 2-8°C.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.