

Goat Anti-Rat IgG heavy and light chain antibody, FITC conjugated

ORDERING INFORMATION

Catalog No.	Size
88204	1.0mg

Format: 1mg/ml in phosphate-buffered saline, 0.2% BSA, 0.09% sodium azide.

SPECIFICATION SUMMARY

Antigen: Rat IgG

Host Species: Goat

Specificity: Antibody is affinity-purified on rat IgG and conjugated to fluorescein isothiocyanate (FITC). Fluorophore: antibody ratio = 5. This antibody may cross-react with IgG of other species.

APPLICATIONS

Immunohistochemistry: use at a dilution of 1:50-1:500.

Immunocytochemistry: use at a dilution of 1:50-1:500.

Immunofluorescence: use at a dilution of 1:50-1:200.

Flow cytometry: use at a dilution of 1:50-1:200.

These are recommended dilutions. Endusers should determine optimal dilutions for their applications. Prepare working dilutions immediately before use.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least 1 year at 2-8°C.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.