

Goat Anti-Rat IgG2c Antibody, FITC Conjugated

ORDERING INFORMATION

| Catalog No. | Size |
|-------------|-------|
| 88503 | 1.0mg |

Format: 1mg/ml in phosphate-buffered saline (PBS), 0.2% BSA, 0.09% sodium azide.

SPECIFICATION SUMMARY

Antigen: Rat IgG2c

Host Species: Goat

Specificity: Antibody is affinity-purified on rat IgG2c and conjugated to fluorescein isothiocyanate. Fluorophore:antibody ratio = 6.2. This antibody may cross-react with IgG2c of other species.

APPLICATIONS

Immunohistochemistry: use at a dilution of 1:50-1:500.

Immunocytochemistry: use at a dilution of 1:50-1:500.

Immunofluorescence: use at a dilution of 1:50-1:500

These are recommended dilutions. Endusers should determine optimal dilutions for their applications. Prepare working dilutions immediately before use.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

This antibody is stable for at least 1 year at 2-8°C.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.